



Using
THE
SAM[®]
SPLINT



The secret is in the **CURVE**



The secret to using the SAM Splint properly is to curve it lengthwise. This longitudinal bend gives the splint its strength.



Curve the outside edges the opposite direction to make it even stronger.



Fold the 9" SAM® Splint around the thumb and over the wrist to stabilize thumb injuries.



Use the SAM® Finger Splint for injured fingers. Fold up the end for extra protection. If you don't have a finger splint, just cut a small part of the SAM® Splint and fold over the exposed edge.



Double the SAM® Splint or create a t-bend for extra strength.



Place the curved splint under the forearm to treat a lower arm fracture or sprain.



Wrap the SAM® Splint around the elbow to splint the wrist or forearm.



For dislocated shoulders, fold the SAM® Splint in threes to form a triangular 'Airplane' splint. Create curves to strengthen.



Using
THE
SAM[®]
JUNIOR



for **ADULT ARMS** and **CHILDREN'S EXTREMITIES**



The 18-inch SAM® JUNIOR will immobilize adult arms and children's extremities. First, curve the splint lengthwise to make it rigid.



For an adult or child's arm, place the SAM® JUNIOR under the forearm to treat a colles fracture or for use as an IV board. Sufferers of carpal tunnel syndrome will find relief wearing this splint at night.



For children's legs, one SAM® Junior can be applied as a sugar-tong splint to immobilize the ankle or leg.



For adults, two SAM® JUNIORs may be used, one applied to each side of the ankle.

Ordering Information

Standard SAM&REG; Splint - 4-1/4" x 36"			
SS-9200-ORG	Orange	\$11.30 ea	\$129.00 cs-12
SS-9200-GRY	Grey		
Junior SAM&REG; Splint - 4-1/4" x 18"			
SS-8200-ORG	Orange	\$7.95 ea	\$174.00 cs-24
SS-8200-GRY	Grey		
Other			
SS-9201	Finger Splint	\$1.00 ea	\$9.50 bg-10
VT-7	Training Video	\$9.95 ea	

© Life-Assist, Inc. All Rights Reserved SAM® Splint is a Registered Trademark TM of The Seaberg Company, Inc.